

# An atomically flat, monocrystalline, gold film thermometer on mica to study energy (heat) exchange at the nano-scale

S. Veronesi<sup>1</sup>, L. Basta<sup>1</sup>, T. Papa<sup>1</sup>, Y. Murata<sup>1</sup>, Z. Dubois<sup>1</sup>, N. Mishra<sup>2,3</sup>, F. Fabbri<sup>2,3</sup>, C. Coletti<sup>2,3</sup>, and S. Heun<sup>1</sup>

<sup>1</sup> NEST, Istituto Nanoscienze-CNR and Scuola Normale Superiore, Piazza S. Silvestro 12, 56127 Pisa, Italy

<sup>2</sup> Center for Nanotechnology Innovation @ NEST, Istituto Italiano di Tecnologia, Piazza S. Silvestro 12, 56127 Pisa, Italy

<sup>3</sup> Graphene Labs, Istituto Italiano di Tecnologia, Via Morego 30, 16163 Genova, Italy

E-mail: stefano.veronesi@nano.cnr.it

There is a great interest in the scientific community to apply calorimetry to samples with mass in the nanogram range. A detailed knowledge of the energy exchange in the fast growing family of micro and nano-systems could allow to obtain valuable information about the chemistry and physics at the nano-scale. The possibility to have an atomically flat thermal probe represents an added value, allowing the unique opportunity to perform Scanning Tunneling Microscopy (STM) experiments together with calorimetry. Here we report on the development of an original experimental setup able to detect temperature variations as low as 10 mK in a sample of  $\sim 10$  ng using a thermometer device having physical dimensions of  $5 \times 5 \text{ mm}^2$ . The technique has been utilized to measure the enthalpy release during the adsorption process of  $\text{D}_2$  on a titanium decorated monolayer graphene (see Fig. 1a). The sensitivity of these thermometers is high enough to detect a hydrogen uptake of  $\sim 10^{-10}$  moles, corresponding to  $\sim 0.2$  ng, with an enthalpy release of about 23  $\mu\text{J}$  [1]. Moreover, we present the fabrication, characterization and calibration of an atomically flat, monocrystalline gold film thermometer on mica substrate. Gold re-crystallization has been obtained inside the STM chamber, allowing the successive investigation of the thermometer surface by LEED and STM imaging (see Fig. 1b). The thermometer calibration demonstrated a thermal exchange coefficient of  $2.2 \times 10^{-7} \text{ W/K}$  and a sensitivity of 1.5 mK at room temperature. The experimental setup allows the simultaneous achievement of heat exchange and surface physics analysis on the same sample.

References:

[1] L. Basta, S. Veronesi, Y. Murata, Z. Dubois, N. P. Mishra, F. Fabbri, C. Coletti and S. Heun, *A sensitive calorimetric technique to study energy (heat) exchange at the nano-scale*, *Nanoscale* **10**, 10079 (2018).




Fig. 1: a) Temperature variation during  $\text{D}_2$  adsorption in Ti functionalized MLG; b) STM image of the Au(111) surface reconstruction.

**Tipo presentazione** (sceglierne uno):                    orale

**Area progettuale di riferimento** (sceglierne una):

SURF Surfaces and interfaces: nanofabrication, imaging, and spectroscopy.